




dualX transforms the fusion environment by delivering a powerful implant through a minimally-invasive approach

Top View



Side View



Description	Est. Graft volume (cm ³)	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Lateral Lumbar Interbody Fusion (LLIF)							
dualX LLIF 0 DEG x 7-9mm x 45mm	3-4 cc	13	22	7	9	54	45
dualX LLIF 0 DEG x 7-9mm x 50mm	4-5 cc	13	22	7	9	59	50
dualX LLIF 0 DEG x 7-9mm x 55mm	5-6 cc	13	22	7	9	64	55
dualX LLIF 7 DEG x 8-11mm x 45mm	3-4 cc	13	22	8	11	54	45
dualX LLIF 7 DEG x 8-11mm x 50mm	4-5 cc	13	22	8	11	59	50
dualX LLIF 7 DEG x 8-11mm x 55mm	5-6 cc	13	22	8	11	64	55
dualX LLIF 7 DEG x 10-13mm x 45mm	4-5 cc	13	22	10	13	54	45
dualX LLIF 7 DEG x 10-13mm x 50mm	5-6 cc	13	22	10	13	59	50
dualX LLIF 7 DEG x 10-13mm x 55mm	6-7 cc	13	22	10	13	64	55
dualX LLIF 7 DEG x 12-15mm x 45mm	6-7 cc	13	22	12	15	54	45
dualX LLIF 7 DEG x 12-15mm x 50mm	7-8 cc	13	22	12	15	59	50
dualX LLIF 7 DEG x 12-15mm x 55mm	8-9 cc	13	22	12	15	64	55

Note: dualX expansion in width and height is definitive and not variable. For example, the dualX LLIF 0 DEG x 7-9mm x 45mm expands exactly from 13mm to 22mm (9mm expansion) in width and 7mm to 9mm (2 mm) in height.



dual expanding interbody fusion system

Top View



Side View



Description	Est. Graft volume (cm ³)	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
Transforaminal Lumbar Interbody Fusion (TLIF)							
dualX TLIF 0 DEG x 7-9mm	2-3 cc	12	21	7	9	39	30
dualX TLIF 4 DEG x 8-10mm	2-3 cc	12	21	8	10	39	30
dualX TLIF 8 DEG x 8-11mm	2-3 cc	12	21	8	11	39	30
dualX TLIF 8 DEG x 9-12mm	2-3 cc	12	21	9	12	39	30
dualX TLIF 8 DEG x 10-13mm	2-3 cc	12	21	10	13	39	30
dualX TLIF 8 DEG x 11-14mm	2-3 cc	12	21	11	14	39	30
dualX TLIF 8 DEG x 12-15mm	2-3 cc	12	21	12	15	39	30

Top View



Side View



Transforaminal/Posterior Lumbar Interbody Fusion (T/PLIF)	Est. Graft volume (cm ³)	Collapsed Width (mm)	Expanded Width (mm)	Collapsed Height (mm)	Expanded Height (mm)	Collapsed Length (mm)	Expanded Length (mm)
dualX T/PLIF 0 DEG x 7-9mm	2-3 cc	12	21	7	9	40	25
dualX T/PLIF 4 DEG x 8-10mm	2-3 cc	12	21	8	10	40	25
dualX T/PLIF 8 DEG x 8-11mm	2-3 cc	12	21	8	11	40	25
dualX T/PLIF 8 DEG x 9-12mm	2-3 cc	12	21	9	12	40	25
dualX T/PLIF 8 DEG x 10-13mm	2-3 cc	12	21	10	13	40	25
dualX T/PLIF 8 DEG x 11-14mm	2-3 cc	12	21	11	14	40	25
dualX T/PLIF 8 DEG x 12-15mm	2-3 cc	12	21	12	15	40	25

Please refer to package insert for current warnings, precautions and instructions for use.
 © Copyright 2019. All Rights Reserved. dualX is a trademark of Amplify Surgical, Inc.
 For more information call 765-AMPLIFY email info@amplifysurgical.com or visit amplifysurgical.com.
 ML10014.A

